

# **VistA Laboratory Enhancements (VLE) – Anatomic Pathology (AP) Order Dialogs**

**Releases: LR\*5.2\*462 and LR\*5.2\*469**

## **Deployment, Installation, Back-Out, and Rollback Guide**



**December 2016**

**Document Version 2.8**

**Department of Veterans Affairs**

**Office of Information and Technology (OI&T)**

## Revision History

Date	Version	Description	Author
12/21/2016	2.8	Minor edits to step 4 of table 5, and figure 1.	REDACTED
12/16/2016	2.7	Minor edits to Appendix A and B.	REDACTED
12/15/2016	2.6	Feedback from Product Support incorporated; re-ordered all of section 4; updated Appendix, edit of figure 2, edit of sections 1.2, 1.3, and 3.	REDACTED
12/14/2016	2.5	Removed IOC Test Site information per feedback from Product Support.	REDACTED
12/12/2016	2.4	Reformatted document.	REDACTED
12/10/2016	2.3	Removed post CPRS version 32 and OR405 instructions.	REDACTED
12/06/2016	2.2	Updated 4.8.2 section for specimen/SNOMED CT mismatch instructions. Removed section 5.2.2. Updated section 3.	REDACTED
11/30/2016	2.1	Addition of new captured installs to Appendix.	REDACTED
11/29/2016	2.0	Removed outdated captured installs.	REDACTED
11/12/2016	1.9	Edit of instructions for Editing Parameter Values; addition of screen captures for section and review provided by Randal Frommater. Appendix A and B updated with captured installs from Test Sites.	REDACTED
11/09/2016	1.8	Addition of instructions for Editing Parameter Values.	REDACTED
09/28/2016	1.7	Updated figure 6 to include SNOMED CT IDs; updated section 4.9.	REDACTED

Date	Version	Description	Author
08/22/2016	1.6	Added notation to Section 4.1.2 step 1c for accession area. Environment check notation added to Section 4.1.2.	REDACTED
07/27/2016	1.5	Removed the Edit Combine Test During Order step as not applicable; removed Figure 2 for How to create a Collection Sample. Added a reference to LEDI IV Installation Guide.	REDACTED
07/25/2016	1.4	Removed the Combine Tests During Ordering Set Step.	REDACTED
07/25/2016	1.3	Addition of Figure 4.	REDACTED
07/20/2016	1.2	Addition of Figure 6; updates to Step 5 in Section 4.1.2.	REDACTED
07/15/2016	1.1	Document baselined.	REDACTED

## Artifact Rationale

This document describes the Deployment, Installation, Back-out, and Rollback Plan for new products going into the VA Enterprise. The plan includes information about system support, issue tracking, escalation processes, and roles and responsibilities involved in all those activities. Its purpose is to provide clients, stakeholders, and support personnel with a smooth transition to the new product or software, and should be structured appropriately, to reflect particulars of these procedures at a single or at multiple locations.

Per the Veteran-focused Integrated Process (VIP) Guide, the Deployment, Installation, Back-out, and Rollback Plan is required to be completed prior to Critical Decision Point #2 (CD #2), with the expectation that it will be updated throughout the lifecycle of the project for each build, as needed.

# Table of Contents

<b>1</b>	<b>INTRODUCTION .....</b>	<b>1</b>
1.1	PURPOSE .....	1
1.2	DEPENDENCIES.....	1
1.3	CONSTRAINTS .....	1
1.4	ORIENTATION.....	2
1.4.1	<i>Computer Dialogue</i> .....	2
1.4.2	<i>Instructions</i> .....	2
1.4.3	<i>User Response</i> .....	2
1.4.4	<i>Warnings and Annotations</i> .....	2
<b>2</b>	<b>ROLES AND RESPONSIBILITIES .....</b>	<b>2</b>
<b>3</b>	<b>DEPLOYMENT.....</b>	<b>3</b>
3.1	TIMELINE .....	4
3.1.1	<i>Site Information (Locations, Deployment Recipients)</i> .....	4
3.1.2	<i>Site Preparation</i> .....	4
3.2	RESOURCES .....	4
3.2.1	<i>Communications</i> .....	4
<b>4</b>	<b>PRE-INSTALLATION AND INSTALLATION PREPARATION INSTRUCTIONS..</b>	<b>4</b>
4.1	PLATFORM INSTALLATION AND PREPARATION.....	5
4.1.1	<i>Installation Timing Recommendation</i> .....	5
4.1.2	<i>Estimated Timing</i> .....	5
4.1.3	<i>Kernel Patches</i> .....	5
4.1.4	<i>Global Growth</i> .....	5
4.2	DOWNLOAD AND EXTRACT FILES .....	5
4.3	DATABASE CREATION .....	5
4.4	INSTALLATION SCRIPTS.....	5
4.5	CRON SCRIPTS .....	6
4.6	ACCESS REQUIREMENTS AND SKILLS NEEDED FOR THE INSTALLATION.....	6
4.7	SYSTEM CONFIGURATION .....	6
4.8	DATABASE TUNING .....	6
4.9	PRE-INSTALLATION AND SYSTEM REQUIREMENTS FOR PATCH LR*5.2*462 .....	6
4.10	INSTALLATION PROCEDURE FOR PATCH LR*5.2*462 .....	6
4.11	PRE-INSTALLATION AND SYSTEM REQUIREMENTS FOR PATCH LR*5.2*469 .....	7
4.12	INSTALLATION PROCEDURE FOR PATCH LR*5.2*469 .....	18
4.13	INSTALLATION VERIFICATION PROCEDURE .....	23
<b>5</b>	<b>BACK-OUT PROCEDURE .....</b>	<b>23</b>
5.1	BACK-OUT STRATEGY .....	23
5.2	BACK-OUT CONSIDERATIONS .....	23
5.2.1	<i>Load Testing</i> .....	23
5.3	BACK-OUT CRITERIA .....	23
5.4	BACK-OUT RISKS.....	24
5.5	AUTHORITY FOR BACK-OUT .....	24

5.6	BACK-OUT PROCEDURE .....	24
5.7	BACK-OUT VERIFICATION PROCEDURE .....	24
<b>6</b>	<b>ROLLBACK PROCEDURE .....</b>	<b>24</b>
6.1	ROLLBACK CONSIDERATIONS .....	24
6.2	ROLLBACK CRITERIA .....	24
6.3	ROLLBACK RISKS .....	24
6.4	AUTHORITY FOR ROLLBACK .....	24
6.5	ROLLBACK PROCEDURE .....	25
6.6	ROLLBACK VERIFICATION PROCEDURE.....	25
<b>APPENDIX A: LR*5.2*462 CAPTURED INSTALLATION .....</b>		<b>26</b>
<b>APPENDIX B: LR*5.2*469 CAPTURED INSTALLATION .....</b>		<b>33</b>

## List of Figures

Figure 1: Examples of AP Tests in the LABORATORY TEST File (#60).....	9
Figure 2: Example of AP Specimen Missing in File #62 .....	11
Figure 3: Example of ORDERABLE ITEMS .....	12
Figure 4: Example of Bone Marrow Display Group Edit.....	15
Figure 5: Example of a Discrepancy in the DISPLAY GROUP for RENAL BIOPSY .....	15
Figure 6: Example of Display Group Edit .....	17
Figure 7: Example of Member Deletion .....	18
Figure 8: Example of Discrepancy for BONE MARROW .....	19
Figure 9: Environment check.....	20
Figure 10: LR*5.2*462 Captured Installation Example .....	26
Figure 11: LR*5.2*469 Captured Installation Example .....	33

## List of Tables

Table 1: Deployment, Installation, Back-out, and Rollback Roles and Responsibilities .....	3
Table 2: Deployment and Installation Schedule .....	4
Table 3: Releases for Installation.....	5
Table 4: M Code Installation Instructions for Patch LR*5.2*462.....	6
Table 5: Pre-Installation Instructions for Patch LR*5.2*469 .....	8
Table 6: M Code Installation Instructions for Patch LR*5.2*469 .....	22

# 1 Introduction

This document describes how to deploy and install patches LR\*5.2\*462 and LR\*5.2\*469, as well as how to back-out the product and rollback to a previous version or data set. This document is a companion to the project charter and management plan for this effort.

Overall, the VLE AP Order Dialog will provide standardization to the current process of AP Order, thereby effectively increasing accuracy, efficiency, and tracking capabilities of specimen processing. The VLE AP Order Dialog will eliminate the manual transcription of handwritten information from the Standard Form (SF) 515 into the existing AP Order Dialog thereby reducing the amount of misidentification of patients, and providing a mechanism to confirm a patient-specimen match.

Patch LR\*5.2\*462 contains a modification to the LABORATORY TEST (#60), LAB ORDER ENTRY (#69), LABORATORY SITE (#69.9) and adds the new file LR CPRS PARAMETERS (#69.71).

Patch LR\*5.2\*469 contains the new file LR AP DIALOG CONFIG (#69.73). The LR AP DIALOG CONFIG (#69.73) file is used to support the building and functioning of the CPRS AP DIALOG for thirteen new tests/orderable items.

The LR AP DIALOG CONFIG file was created to support the new LR OTHER LAB AP TESTS order dialog.

## 1.1 Purpose

The purpose of this guide is to provide a single, common document that describes how, when, where, and to whom patches LR\*5.2\*462 and LR\*5.2\*469 will be deployed and installed, as well as how it is to be backed out and rolled back, if necessary. The plan also identifies resources, communications plan, and rollout schedule. Specific instructions for installation, back-out, and rollback are included in this document.

## 1.2 Dependencies

Patch LR\*5.2\*462 contains the following dependencies:

- LR\*5.2\*259
- LR\*5.2\*291
- LR\*5.2\*350
- LR\*5.2\*422
- LR\*5.2\*433
- LR\*5.2\*450

Patch LR\*5.2\*469 contains the following dependency:

- LR\*5.2\*462.

## 1.3 Constraints

Security controls will be inherited from VistA and therefore will be fully compliant with National Institute of Standards and Technology (NIST) controls and in compliance with

Directive 6500. In addition, patches LR\*5.2\*462 and LR\*5.2\*469 will be 508 compliant and designed to ensure no performance impacts will be experienced in the production environments.

## 1.4 Orientation

This section addresses package or audience specific notations or directions (e.g., symbols used to indicate terminal dialogues or user responses) for the installation and post-installation instructions included in this document.

All headings and text in this guide are intentionally formatted flush left, regardless of the heading level, to save space and to make for better readability.

In tables which list mandatory steps (for pre-installation and installation), a column is provided at the right-hand side so that users may check (✓) off the step as it is performed.

### 1.4.1 Computer Dialogue

The computer dialogue will appear in Courier New 11-point font.

**Example:** `Courier New font 11 points`

### 1.4.2 Instructions

Instructions will appear in Arial 11-point font.

**Example:** Arial font 11 points

### 1.4.3 User Response


User entry responses will appear in Courier New 11-point font.

**Example:** `Courier New font 11 points`

### 1.4.4 Warnings and Annotations

The verbiage warning and note will appear in Arial 11-point bold font to alert the user about important information that will follow.

**Example: Arial bold font 11 points**

Warnings of high importance will be denoted with a  symbol.

## 2 Roles and Responsibilities

The deployment, installation, back-out, and roll-back roles and responsibilities are outlined in table 1 which also depicts the teams that will perform the steps described in this plan.

**Table 1: Deployment, Installation, Back-out, and Rollback Roles and Responsibilities**

ID	Team	Phase / Role	Tasks
1	Field Operations (FO)	Deployment	Plan and schedule deployment (including orchestration with vendors).
2	FO	Deployment	Determine and document the roles and responsibilities of those involved in the deployment.
3	FO	Deployment	Test for operational readiness.
4	FO	Deployment	Execute deployment.
5	FO	Installation	Plan and schedule installation.
6	FO	Installation	Ensure authority to operate and that certificate authority security documentation is in place.
7	FO	Installation	Validate through facility POC to ensure that IT equipment has been accepted using asset inventory processes.
8	FO	Installations	Coordinate training.
9	FO	Back-out	Confirm availability of back-out instructions and back-out strategy (what are the criteria that trigger a back-out).
10	FO	Post Deployment	Hardware, Software and System Support.

### 3 Deployment

A National Release is planned for December 2016 for LR\*5.2\*462 and LR\*5.2\*469.

Software Development support will be available for forty-five days after the software has been released nationally.

Deployment will be performed by Local Facility staff and supported by team members from one or more of the operations organizations: Enterprise Systems Engineering (ESE), Field Operations (FO), Enterprise Operations (EO), and/or others.



## 3.1 Timeline

A high level overview of the Installation schedule with estimated dates is shown in table 2.

**Table 2: Deployment and Installation Schedule**

Build	Site	Date of Installation	Target
National Release	National	December 2016	Production at all sites for patches LR*5.2*462 and LR*5.2*469.

### 3.1.1 Site Information (Locations, Deployment Recipients)

Patches LR\*5.2\*462 and LR\*5.2\*469 were tested by the following sites:

- REDACTED
- REDACTED
- REDACTED
- REDACTED
- REDACTED

### 3.1.2 Site Preparation

Patches LR\*5.2\*462 and LR\*5.2\*469 are enhancements to the existing VistA and CPRS legacy applications; therefore, new equipment is not required.

## 3.2 Resources

LR\*5.2\*462 and LR\*5.2\*469 are enhancements to the existing VistA and CPRS legacy applications; therefore, new equipment is not required.

### 3.2.1 Communications

Technicians will use email and conference calls for communication during the deployment; Email and/or conference calls will be utilized to notify the stakeholders of the successful release of the product.

## 4 Pre-Installation and Installation Preparation Instructions

Pre-installation and installation instructions for patches LR\*5.2\*462 and LR\*5.2\*469 are provided sequentially. Follow the installation instructions for patch LR\*5.2\*462 first followed by the instructions for the installation of patch LR\*5.2\*469 second.

## 4.1 Platform Installation and Preparation

The Patch installation needs to be coordinated with the LIM at each site.

### 4.1.1 Installation Timing Recommendation

This patch may be loaded with users on the system. If desired, perform the installation during non-peak hours.

### 4.1.2 Estimated Timing

The installation will take less than 5 minutes.

### 4.1.3 Kernel Patches

Kernel patches must be current on the target system to avoid problems loading and/or installing this patch.

### 4.1.4 Global Growth

There is no significant change to global growth.

## 4.2 Download and Extract Files

Installation of the builds shall be achieved via Packman messages. Associated documentation will be available on the VistA Document Library (VDL). The online versions will be updated as needed. Please look for the latest version on the VDL:

<http://www.va.gov/vdl>

The following table lists the releases that are available for installation.

**Table 3: Releases for Installation**

Host File	File Name
Packman Distribution	Each Packman message contains 1 patch: <ul style="list-style-type: none"><li>• LR*5.2*462</li><li>• LR*5.2*469</li></ul>

## 4.3 Database Creation

This section is not applicable as there will not be any databases created.

## 4.4 Installation Scripts

This section is not applicable as there will not be any installation scripts.

## 4.5 CRON Scripts

This section is not applicable as there will not be any CRON scripts.

## 4.6 Access Requirements and Skills Needed for the Installation

This guide is written with the assumption that the reader is experienced or familiar with the following:

- VistA computing environment
- VA FileMan data structures and terminology
- Microsoft Windows

Product Support will coordinate the patch installation with the LIM at each site. The LIM and Product Support have the necessary access and skill set to conduct the installation.

## 4.7 System Configuration

The pre-installations instructions complete the system configuration process.

## 4.8 Database Tuning

Database tuning is not applicable.

## 4.9 Pre-installation and System Requirements for Patch LR\*5.2\*462

Patch LR\*5.2\*462 will install new fields to the LABORATORY TEST (#60), LAB ORDER ENTRY (#69) and LABORATORY SITE PARAMETERS (#69.9) files. It will also install a new file, LR CPRS SCREEN (#69.71), four new "LR7OAP" options, a print template for the LR CPRS SCREEN file and an edit template for the LABORATORY TEST file.

## 4.10 Installation Procedure for Patch LR\*5.2\*462

Follow the instructions outlined in the table below to perform the installation of Patch LR\*5.2\*462.

**Table 4: M Code Installation Instructions for Patch LR\*5.2\*462**

Step #	Description	✓
1	On the PackMan menu, use the <code>INSTALL/CHECK MESSAGE</code> option. This option loads the patch into a Transport Global on your system.	
2	From the Kernel Installation and Distribution System (KIDS) menu, select the <code>Installation</code> menu.	
3	From this menu, you may elect to use the following options:	

VLE AP LR\*5.2\*462, LR\*5.2\*469

Step #	Description		✓
	<b>a</b>	Verify Checksums in Transport Global: This option will allow you to ensure the integrity of the routines that are in the transport global.	
	<b>b</b>	Print Transport Global: This option will allow you to view the components of the KIDS build.	
	<b>c</b>	Compare Transport Global to Current System: This option will allow you to view all changes that will be made when this patch is installed. It compares all components of this patch (routines, Data Dictionaries (DD's), templates, etc.).	
	<b>d</b>	Backup a Transport Global: This option will create a backup message of any routines exported with this patch. It will not backup any other changes such as DD's or templates.	
<b>4</b>	Under the Installation menu, select the package LR*5.2*462.		
<b>5</b>	When prompted 'Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//', respond "NO".		
<b>6</b>	When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//', respond "NO".		
<b>7</b>	When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//', respond "NO".		
<b>8</b>	Enter the Device you want to print the Install message. You can queue the install by enter a "Q" at the device prompt. Enter a "^" to abort the install.		
<b>9</b>	When prompted 'Delay Install (Minutes): (0 - 60) ://', respond "0".		

## 4.11 Pre-installation and System Requirements for Patch LR\*5.2\*469

After performing the pre-installation instructions and the installation of patch LR\*5.2\*462, perform the pre-installation instructions and installation of patch LR\*5.2\*469.

The LIM must create thirteen new lab tests to support the CPRS AP DIALOG functionality. Before the installation of the patch, an environment check will take place to ensure that the following criteria have been met: the orderable item name has been properly identified with a Display Group of Anatomic Pathology, and the preferred SNOMED CT code for the specimens have been properly identified. If all of the criteria have been met, the installation will proceed; otherwise, the installation will fail and the LIM will need to rectify any detected discrepancy.



**Note:** In the event that a discrepancy occurs, due to a mismatch for the specimen name or the SNOMED CT, contact the National Service Desk (NSD) at 855-673-4357 to log an incident ticket.

Follow the instructions outlined in the following table for Patch LR\*5.2\*469.

**Table 5: Pre-Installation Instructions for Patch LR\*5.2\*469**

Step #	Description	✓
1	The Automated Data Processing Application Coordinator (ADPAC) at your facility will need to complete this step.	
	<p><b>a</b></p> <p>Create thirteen new lab tests with a collection sample of AP SPECIMEN by using the appropriate institution entries for your facility for the following lab tests/orderable items:</p> <ul style="list-style-type: none"><li>• BONE MARROW</li><li>• BRONCHIAL BIOPSY</li><li>• BRONCHIAL CYTOLOGY</li><li>• DERMATOLOGY</li><li>• FINE NEEDLE ASPIRATE</li><li>• GASTROINTESTINAL ENDOSCOPY</li><li>• GENERAL FLUID</li><li>• GYNECOLOGY (PAP SMEAR)</li><li>• RENAL BIOPSY</li><li>• TISSUE EXAM</li><li>• URINE</li><li>• UROLOGY, PROSTATE</li><li>• UROLOGY, BLADDER/URETER</li></ul> <p><b>WARNING:</b> Ensure that there are no other entries in the LABORATORY TEST file that have the same lab test name or synonym even if the subscript is different.</p>	

Step #	Description	✓
1	Step 1 Continued...	
	<b>b</b> Set up all tests with a Type of Neither. <b>Note:</b> See figure 1 for examples of how the lab tests may appear in the LABORATORY TEST file (#60) after all of the entries have been created. <b>Note:</b> Workload code configuration should be done in compliance with the facility's current laboratory practices.	
	<b>c</b> <b>WARNING:</b> If the COLLECTION SAMPLE file (#62) does not have an AP SPECIMEN entry, follow the instructions outlined in Section 2.2 of the LEDI IV Installation Guide or call Product Support for assistance. An example of a warning dialog for the AP Specimen entry is shown in figure 2 below. <b>Note:</b> There is no restriction on what is utilized as the Accession Area.	

**Figure 1: Examples of AP Tests in the LABORATORY TEST File (#60)**

NAME: BONE MARROW	TYPE: NEITHER
SUBSCRIPT: SURGICAL PATHOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: BMARROW	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: SURGICAL PATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Surgical Pathology Tissue Exam	
CPRS SCREEN: 118-BONE MARROW	
NAME: BRONCHIAL BIOPSY	TYPE: NEITHER
SUBSCRIPT: SURGICAL PATHOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: BROBIOSP	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: SURGICAL PATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Surgical Pathology Tissue Exam	
CPRS SCREEN: 109-BRONCHIAL BIOPSY	
NAME: BRONCHIAL CYTOLOGY	TYPE: NEITHER
SUBSCRIPT: CYTOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: BRONC-CY	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: CYTOPATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Report Cytology	
CPRS SCREEN: 102-BRONCHIAL CYTOLOGY	

NAME: GYNECOLOGY (PAP SMEAR) TYPE: NEITHER

SUBSCRIPT: CYTOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: GYN-PAP	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: CYTOPATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Report Cytology	
CPRS SCREEN: 126-GYNECOLOGY (PAP SMEAR)	

NAME: URINE	TYPE: NEITHER
SUBSCRIPT: CYTOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: URCYTO	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: CYTOPATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Report Cytology	
CPRS SCREEN: 127-URINE	

NAME: UROLOGY,BLADDER/URETER	TYPE: NEITHER
SUBSCRIPT: SURGICAL PATHOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: UROBLURE	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: SURGICAL PATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Surgical Pathology Tissue Exam	
CPRS SCREEN: 123-UROLOGY,BLADDER/URETER	

NAME: UROLOGY,PROSTATE	TYPE: NEITHER
SUBSCRIPT: SURGICAL PATHOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: UROPROST	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: SURGICAL PATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Surgical Pathology Tissue Exam	
CPRS SCREEN: 121-UROLOGY,PROSTATE	

NAME: RENAL BIOPSY	TYPE: NEITHER
SUBSCRIPT: SURGICAL PATHOLOGY	HIGHEST URGENCY ALLOWED: STAT
REQUIRED COMMENT: ORDER COMMENT	COMBINE TEST DURING ORDER: NO
PRINT NAME: RENAL BX	
COLLECTION SAMPLE: AP SPECIMEN	
INSTITUTION: CAMP MASTER	ACCESSION AREA: SURGICAL PATHOLOGY
SITE NOTES DATE: MAY 03, 2016	
NOTE: Test added in support of the VLE-AP order dialog enhancement	
NATIONAL VA LAB CODE: Surgical Pathology Tissue Exam	
CPRS SCREEN: 128-RENAL BIOPSY	

**Figure 2: Example of AP Specimen Missing in File #62**

\*\*\* 'ABDOMEN'



```

Found SPECIMEN for SNOMED CT '113345001'

Additional Specimens found but not listed in this example.

***** AP SPECIMEN *****
* COLLECTION SAMPLE (File 62) 'AP SPECIMEN' is missing! *
***** AP SPECIMEN *****

LR*5.2*469 Build will not be installed, Transport Global deleted!
Dec 07, 2016@11:01:26

```

Step #	Description	✓
2	Confirm that the lab tests were created properly in the ORDERABLE ITEMS file (#101.43). The number seen in the ID field will point to the internal entry number (IEN) for the test in your LABORATORY TEST file (#60). See figure 3 for an example of how the ORDERABLE ITEMS might appear in the ORDERABLE ITEMS file (#101.43).	

**Figure 3: Example of ORDERABLE ITEMS**

NAME: BONE MARROW	ID: 719;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: BONE MARROW
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: BRONCHIAL BIOPSY	ID: 5125;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: BRONCHIAL BIOPSY
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: BRONCHIAL CYTOLOGY	ID: 5126;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: BRONCHIAL CYTOLOGY
SET: LAB	
SET: AP	
LAB SECTION: CYTOLOGY	TYPE: neither
NAME: DERMATOLOGY	ID: 5127;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: DERMATOLOGY
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: FINE NEEDLE ASPIRATE	ID: 5128;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: FINE NEEDLE ASPIRATE
SET: LAB	
SET: AP	
LAB SECTION: CYTOLOGY	TYPE: neither

NAME: GENERAL FLUID	ID: 5131;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: GENERAL FLUID
SET: LAB	
SET: AP	
LAB SECTION: CYTOLOGY	TYPE: neither
NAME: TISSUE EXAM	ID: 5133;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: TISSUE EXAM
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: GASTROINTESTINAL ENDOSCOPY	ID: 5129;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: GASTROINTESTINAL
ENDOSCOPY	
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: GYNECOLOGY (PAP SMEAR)	ID: 5132;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: GYNECOLOGY (PAP
SMEAR)	
SET: LAB	
SET: AP	
LAB SECTION: CYTOLOGY	TYPE: neither
NAME: URINE	ID: 5137;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: URINE
SET: LAB	
SET: AP	
LAB SECTION: CYTOLOGY	TYPE: neither
NAME: UROLOGY,BLADDER/URETER	ID: 5134;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME:
UROLOGY,BLADDER/URETER	
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: UROLOGY,PROSTATE	ID: 5135;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: UROLOGY,PROSTATE
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither
NAME: RENAL BIOPSY	ID: 5136;99LRT
DISPLAY GROUP: LABORATORY	PACKAGE NAME: RENAL BIOPSY
SET: LAB	
SET: AP	
LAB SECTION: SURGICAL PATHOLOGY	TYPE: neither

Step #	Description	✓
3	<p>Edit the ORDERABLE ITEM file (#101.43) and change the DISPLAY GROUP from LABORATORY to ANATOMIC PATHOLOGY for each of the following thirteen lab tests:</p> <ul style="list-style-type: none"> <li>• BONE MARROW</li> <li>• BRONCHIAL BIOPSY</li> <li>• BRONCHIAL CYTOLOGY</li> <li>• DERMATOLOGY</li> <li>• FINE NEEDLE ASPIRATE</li> <li>• GASTROINTESTINAL ENDOSCOPY</li> <li>• GENERAL FLUID</li> <li>• GYNECOLOGY (PAP SMEAR)</li> <li>• RENAL BIOPSY</li> <li>• TISSUE EXAM</li> <li>• URINE</li> <li>• UROLOGY, PROSTATE</li> <li>• UROLOGY, BLADDER/URETER</li> </ul> <p><b>Note:</b> This edit will remove the LAB entry from the SET field.  <b>Note:</b> See figure 4 below for an example.</p>	
	<p><b>a</b> When prompted, 'Select OPTION:' respond 'ENTER OR EDIT FILE ENTRIES'</p>	
	<p><b>b</b> When prompted, 'INPUT TO WHAT FILE:' respond 'ORDERABLE ITEMS'</p>	
	<p><b>c</b> When prompted, 'EDIT WHICH FIELD:' respond 'ALL// DISPLAY GROUP'</p>	
	<p><b>d</b> When prompted 'Select ORDERABLE ITEMS NAME:' respond with the name of each of the thirteen lab tests.</p>	
	<p><b>e</b> When prompted 'DISPLAY GROUP:' respond 'LABORATORY// ANATOMIC PATHOLOGY'</p>	
	<p><b>f</b> When prompted 'Select ORDERABLE ITEMS NAME:' respond with the name of each of the thirteen lab tests.</p> <p><b>Note:</b> You must change the DISPLAY GROUP from LABORATORY to ANATOMIC PATHOLOGY for each of the thirteen lab tests. See figure 4 below for an example.</p>	

**Figure 4: Example of Bone Marrow Display Group Edit**

```
Select OPTION: ENTER OR EDIT FILE ENTRIES

INPUT TO WHAT FILE: ORDERABLE ITEMS//
EDIT WHICH FIELD: ALL// DISPLAY GROUP
THEN EDIT FIELD:

Select ORDERABLE ITEMS NAME:      BONE MARROW
DISPLAY GROUP: LABORATORY// ANATOMIC PATHOLOGY

Select ORDERABLE ITEMS NAME:

Review of Orderable Item:
  NAME: BONE MARROW                      ID: 719;99LRT
  DISPLAY GROUP: ANATOMIC PATHOLOGY      PACKAGE NAME: BONE MARROW
  SET: AP
  LAB SECTION: SURGICAL PATHOLOGY        TYPE: neither
```

The following figure is an example of the DISPLAY GROUP that has not been properly changed which was identified during the environment check. The user will need to correct the discrepancy by following step #3 outlined above.

**Figure 5: Example of a Discrepancy in the DISPLAY GROUP for RENAL BIOPSY**

```
Found ORDERABLE ITEM 'BONE MARROW'.
Found ORDERABLE ITEM 'BRONCHIAL BIOPSY'.
Found ORDERABLE ITEM 'BRONCHIAL CYTOLOGY'.
Found ORDERABLE ITEM 'DERMATOLOGY'.
Found ORDERABLE ITEM 'FINE NEEDLE ASPIRATE'.
Found ORDERABLE ITEM 'GASTROINTESTINAL ENDOSCOPY'.
Found ORDERABLE ITEM 'GENERAL FLUID'.
Found ORDERABLE ITEM 'GYNECOLOGY (PAP SMEAR)'.
Found ORDERABLE ITEM 'TISSUE EXAM'.
Found ORDERABLE ITEM 'UROLOGY,BLADDER/URETER'.
Found ORDERABLE ITEM 'UROLOGY,PROSTATE'.
*** There is no ORDERABLE ITEM 'RENAL BIOPSY' within the ANATOMIC PATHOLOGY
display group

LR*5.2*469 Build will not be installed, Transport Global deleted!
Current date and time will be displayed here.
```

Step #	Description	✓
4	<p><b>Note:</b> This step will need to be coordinated with your Clinical Coordinator and/or IT Support depending on access rights.</p> <p>Edit the DISPLAY GROUPS file (#100.98) and add SURGICAL PATHOLOGY, CYTOLOGY, and ELECTRON MICROSCOPY as new members. This edit will display the anatomic pathology orders under the ANATOMIC PATHOLOGY DISPLAY GROUP instead of the LABORATORY DISPLAY GROUP within the CPRS Order Tab.</p>	
	<b>a</b> When prompted 'INPUT TO WHAT FILE:' respond 'DISPLAY GROUP'	
	<b>b</b> When prompted 'EDIT WHICH FIELD:' respond 'ALL'	
	<b>c</b> When prompted 'Select DISPLAY GROUP NAME:' respond 'ANATOMIC PATHOLOGY'	
	<b>d</b> When prompted 'NAME:' respond 'ANATOMIC PATHOLOGY'	
	<b>e</b> When prompted 'Select MEMBER:' respond 'SURGICAL PATHOLOGY'	
	<b>f</b> When prompted 'Are you adding 'SURGICAL PATHOLOGY' as a new MEMBER (the 1ST for this DISPLAY GROUP)?' respond 'YES'	
	<b>g</b> When prompted 'MEMBER SEQUENCE: 1//' respond with the ENTER key.	
	<b>h</b> When prompted 'Select MEMBER:' respond 'CYTOLOGY'	
	<b>i</b> When prompted 'Are you adding 'CYTOLOGY' as a new MEMBER (the 2ND for this DISPLAY GROUP)?' respond 'YES'	
	<b>j</b> When prompted 'MEMBER SEQUENCE: 2//' respond with the ENTER key.	
	<b>k</b> When prompted 'Select MEMBER:' respond 'ELECTRON MICROSCOPY'	
	<b>l</b> When prompted 'Are you adding 'ELECTRON MICROSCOPY' as a new MEMBER (the 3RD for this DISPLAY GROUP)?' respond 'YES'	
	<b>m</b> When prompted 'MEMBER SEQUENCE: 3//' respond with the ENTER key.	
	<b>n</b> When prompted 'MIXED NAME: Anat. Path.//' respond with the ENTER key.	
	<b>o</b> When prompted 'SHORT NAME: AP//' respond with the ENTER key.	
	<b>p</b> When prompted 'DEFAULT DIALOG' respond with the ENTER key to accept the default of a blank entry.	

**Figure 6: Example of Display Group Edit**

```

INPUT TO WHAT FILE: DISPLAY GROUP//
EDIT WHICH FIELD: ALL//

Select DISPLAY GROUP NAME: ANATOMIC PATHOLOGY
NAME: ANATOMIC PATHOLOGY//
Select MEMBER: SURGICAL PATHOLOGY
Are you adding 'SURGICAL PATHOLOGY' as a new MEMBER (the 1ST for
this DISPLAY GROUP)? No// Y (Yes)
MEMBER SEQUENCE: 1//
Select MEMBER: CYTOLOGY
Are you adding 'CYTOLOGY' as a new MEMBER (the 2ND for this
DISPLAY GROUP)? No // Y (Yes)
MEMBER SEQUENCE: 2//
Select MEMBER: ELECTRON MICROSCOPY
Are you adding 'ELECTRON MICROSCOPY' as a new MEMBER (the 3RD for
this DISPLAY GROUP)? No// Y (Yes)
MEMBER SEQUENCE: 3//
Select MEMBER:
MIXED NAME: Anat. Path.//
SHORT NAME: AP//
DEFAULT DIALOG:[leave blank]

```

Step #	Description	✓
5	<p><b>Note:</b> This step will need to be coordinated with your Clinical Coordinator.</p> <p>Delete the AP display groups SURGICAL PATHOLOGY, CYTOLOGY, and ELECTRON MICROSCOPY from the LABORATORY display group entry.</p>	
	<b>a</b> When prompted 'Select DISPLAY GROUP NAME:' respond 'LABORATORY'	
	<b>b</b> When prompted 'Select MEMBER:' respond 'SURGICAL PATHOLOGY'. To delete the member, enter the '@' sign.	
	<b>c</b> When prompted 'SURE YOU WANT TO DELETE THE ENTIRE MEMBER?' respond 'Yes'	
	<b>d</b> When prompted 'Select MEMBER:' respond 'CYTOLOGY'. To delete the member, enter the '@' sign.	
	<b>e</b> When prompted 'SURE YOU WANT TO DELETE THE ENTIRE MEMBER?' respond 'Yes'	
	<b>f</b> When prompted 'Select MEMBER:' respond 'ELECTRON MICROSCOPY'. To delete the member, enter the '@' sign.	
	<b>g</b> When prompted 'SURE YOU WANT TO DELETE THE ENTIRE MEMBER?' respond 'Yes'	

**Figure 7: Example of Member Deletion**

```
Select DISPLAY GROUP NAME: lab  LABORATORY
NAME: LABORATORY//
Select MEMBER: CYTOLOGY// @
    SURE YOU WANT TO DELETE THE ENTIRE MEMBER? Y  (Yes)
Select MEMBER: AUTOPSY// SURGICAL PATHOLOGY
    ...OK? Yes//  (Yes)

    MEMBER: SURGICAL PATHOLOGY// @
    SURE YOU WANT TO DELETE THE ENTIRE MEMBER? Y  (Yes)
Select MEMBER: ELECT  RON MICROSCOPY
    ...OK? Yes//  (Yes)

    MEMBER: ELECTRON MICROSCOPY// @
    SURE YOU WANT TO DELETE THE ENTIRE MEMBER? Y  (Yes)
Select MEMBER: ?
    Answer with MEMBER, or SEQUENCE
    Choose from:
    1          CHEMISTRY
    2          HEMATOLOGY
    3          MICROBIOLOGY
    4          BLOOD BANK
    5          ANATOMIC PATHOLOGY
    8          AUTOPSY

    You may enter a new MEMBER, if you wish
    Enter a display group.  A display group that is an ancestor
    may not also be a member.
    Answer with DISPLAY GROUP NAME
    Do you want the entire 59-Entry DISPLAY GROUP List?

    You may enter a new DISPLAY GROUP, if you wish
    ANSWER MUST BE 3-30 CHARACTERS IN LENGTH
Select MEMBER:
MIXED NAME: Lab//
SHORT NAME: LAB//
DEFAULT DIALOG: LR OTHER LAB TESTS//
```

## 4.12 Installation Procedure for Patch LR\*5.2\*469

Before the installation of patch LR\*5.2\*469 takes place, an environment check will be performed to ensure that the ORDERABLE ITEMS (from File #101.43) and their associated TOPOGRAPHY/SNOMED CT ID are an exact match and that multiple instances do not exist. If an exact match is not found or if there are multiple instances, the user will be informed and the environmental check will fail; the user will need to correct the discrepancy before performing the installation. The following figure shows an example of a discrepancy for BONE MARROW.



**Note:** If the Environment Check fails due to a mismatch for the specimen name or the SNOMED CT, contact the NSD at 855-673-4357 to log an incident ticket.

**Figure 8: Example of Discrepancy for BONE MARROW**

***	Unable to find the ORDERABLE ITEM 'BONE MARROW'
Found	ORDERABLE ITEM 'BRONCHIAL CYTOLOGY'
Found	ORDERABLE ITEM 'BRONCHIAL BIOPSY'
Found	ORDERABLE ITEM 'DERMATOLOGY'
Found	ORDERABLE ITEM 'FINE NEEDLE ASPIRATE'
Found	ORDERABLE ITEM 'GENERAL FLUID'
Found	ORDERABLE ITEM 'TISSUE EXAM'
Found	ORDERABLE ITEM 'GASTROINTESTINAL ENDOSCOPY'
Found	ORDERABLE ITEM 'GYNECOLOGY (PAP SMEAR)'
Found	ORDERABLE ITEM 'URINE'
Found	ORDERABLE ITEM 'UROLOGY, PROSTATE'
Found	ORDERABLE ITEM 'UROLOGY, BLADDER/URETER'
Found	ORDERABLE ITEM 'RENAL BIOPSY'
Found	SPECIMEN 'ABDOMEN'
Found	SPECIMEN 'ADRENAL GLAND'
Found	SPECIMEN 'ANAL CANAL'
Found	SPECIMEN 'APPENDIX'
Found	SPECIMEN 'ASCENDING COLON'
Found	SPECIMEN 'BILE DUCT CYTOLOGIC MATERIAL'
Found	SPECIMEN 'BILE DUCT MUCOUS MEMBRANE'
Found	SPECIMEN 'BILIARY TRACT'
Found	SPECIMEN 'BODY OF PANCREAS'
Found	SPECIMEN 'BONE'
***	Unable to find the SPECIMEN 'BONE MARROW'
Found	SPECIMEN 'BREAST'
Found	SPECIMEN 'BRONCHIAL CYTOLOGIC MATERIAL'
Found	SPECIMEN 'BRONCHUS'
Found	SPECIMEN 'CARDIAC INCISURE OF STOMACH'
Found	SPECIMEN 'CARDIAC OSTIUM OF STOMACH'
Found	SPECIMEN 'CARDIO-ESOPHAGEAL JUNCTION'
Found	SPECIMEN 'CECUM'
Found	SPECIMEN 'CEREBROSPINAL FLUID CYTOLOGIC MATERIAL'
Found	SPECIMEN 'CERVICAL CYTOLOGIC MATERIAL'
Found	SPECIMEN 'CHORIONIC VILLI'
Found	SPECIMEN 'COLON'
Found	SPECIMEN 'DESCENDING COLON'
Found	SPECIMEN 'DUODENUM'
Found	SPECIMEN 'ESOPHAGUS'
Found	SPECIMEN 'ESOPHAGUS, LOWER THIRD'
Found	SPECIMEN 'ESOPHAGUS, MIDDLE THIRD'
Found	SPECIMEN 'ESOPHAGUS, UPPER THIRD'
Found	SPECIMEN 'GASTRIC FUNDUS'
Found	SPECIMEN 'GASTRIC JUICE'
Found	SPECIMEN 'GREATER CURVATURE OF STOMACH'
Found	SPECIMEN 'HEAD OF PANCREAS'
Found	SPECIMEN 'ILEUM'
Found	SPECIMEN 'JEJUNUM'
Found	SPECIMEN 'KIDNEY'
Found	SPECIMEN 'LEFT COLIC FLEXURE'
Found	SPECIMEN 'LEFT TESTIS'
Found	SPECIMEN 'LESSER CURVATURE OF STOMACH'
Found	SPECIMEN 'LIVER'
Found	SPECIMEN 'LUNG'
Found	SPECIMEN 'LYMPH NODE'
Found	SPECIMEN 'MEDIASTINUM'
Found	SPECIMEN 'NECK, LEFT SIDE'
Found	SPECIMEN 'NECK, RIGHT SIDE'
Found	SPECIMEN 'PANCREAS'
Found	SPECIMEN 'PAROTID GLAND'
Found	SPECIMEN 'PELVIS'
Found	SPECIMEN 'PERICARDIAL CYTOLOGIC MATERIAL'



```

Found SPECIMEN 'PERIRENAL TISSUE'
Found SPECIMEN 'PERITONEAL CYTOLOGIC MATERIAL'
Found SPECIMEN 'PERITONEUM'
Found SPECIMEN 'PLEURAL CYTOLOGIC MATERIAL'
Found SPECIMEN 'PROSTATE'
Found SPECIMEN 'PYLORIC ANTRUM'
Found SPECIMEN 'PYLORUS'
Found SPECIMEN 'RECTOSIGMOID JUNCTION'
Found SPECIMEN 'RECTUM'
Found SPECIMEN 'RENAL PELVIS'
Found SPECIMEN 'RETROPERITONEUM'
Found SPECIMEN 'RIGHT COLIC FLEXURE'
Found SPECIMEN 'RIGHT TESTIS'
Found SPECIMEN 'SALIVARY GLAND'
Found SPECIMEN 'SIGMOID COLON'
Found SPECIMEN 'SKIN'
Found SPECIMEN 'SMALL INTESTINE'
Found SPECIMEN 'SOFT TISSUES'
Found SPECIMEN 'SPLEEN'
Found SPECIMEN 'STOMACH'
Found SPECIMEN 'SUBLINGUAL GLAND'
Found SPECIMEN 'SUBMAXILLARY GLAND'
Found SPECIMEN 'SUBPHRENIC FOSSA'
Found SPECIMEN 'SYNOVIAL CYTOLOGIC MATERIAL'
Found SPECIMEN 'TAIL OF PANCREAS'
Found SPECIMEN 'THYROID GLAND'
Found SPECIMEN 'TISSUE'
Found SPECIMEN 'TRANSVERSE COLON'
Found SPECIMEN 'URETER'
Found SPECIMEN 'URINARY BLADDER'
Found SPECIMEN 'URINE'
Found SPECIMEN 'VAS DEFERENS'
Found SPECIMEN 'VERTEBRA'

```

The following figure shows what is verified during the environment check for the ORDERABLE ITEMS and associated TOPOGRAPHY/SNOMED CT ID.



**Note:** If the Environment Check fails due to a mismatch for the specimen name or the SNOMED CT, contact the NSD at 855-673-4357 to log an incident ticket.

**Figure 9: Environment check**

ORDERABLE ITEMS checked:

```

BONE MARROW
BRONCHIAL CYTOLOGY
BRONCHIAL BIOPSY
DERMATOLOGY
FINE NEEDLE ASPIRATE
GENERAL FLUID
TISSUE EXAM
GASTROINTESTINAL ENDOSCOPY
GYNECOLOGY (PAP SMEAR)
URINE
UROLOGY, PROSTATE
UROLOGY, BLADDER/URETER
RENAL BIOPSY

```

TOPOGRAPHY(specimens)/SNOMED CT ID codes checked:

ABDOMEN - 113345001  
 ADRENAL GLAND - 23451007  
 ANAL CANAL - 34381000  
 APPENDIX - 66754008  
 ASCENDING COLON - 9040008  
 BILE DUCT CYTOLOGIC MATERIAL - 110928002  
 BILE DUCT MUCOUS MEMBRANE - 7035006  
 BILIARY TRACT - 34707002  
 BODY OF PANCREAS - 40133006  
 BONE - 3138006  
 BONE MARROW - 14016003  
 BREAST - 76752008  
 BRONCHIAL CYTOLOGIC MATERIAL - 110912007  
 BRONCHUS - 955009  
 CARDIAC INCISURE OF STOMACH - 5459006  
 CARDIAC OSTIUM OF STOMACH - 63853002  
 CARDIO-ESOPHAGEAL JUNCTION - 25271004  
 CECUM - 32713005  
 CEREBROSPINAL FLUID CYTOLOGIC MATERIAL - 110969006  
 CERVICAL CYTOLOGIC MATERIAL - 110949001  
 CHORIONIC VILLI - 2049008  
 COLON - 71854001  
 DESCENDING COLON - 32622004  
 DUODENUM - 38848004  
 ESOPHAGUS - 32849002  
 ESOPHAGUS, LOWER THIRD - 67173009  
 ESOPHAGUS, MIDDLE THIRD - 19000002  
 ESOPHAGUS, UPPER THIRD - 54738009  
 GASTRIC FUNDUS - 414003  
 GASTRIC JUICE - 31773000  
 GREATER CURVATURE OF STOMACH - 89382009  
 HEAD OF PANCREAS - 64163001  
 ILEUM - 34516001  
 JEJUNUM - 21306003  
 KIDNEY - 64033007  
 LEFT COLIC FLEXURE - 72592005  
 LEFT TESTIS - 63239009  
 LESSER CURVATURE OF STOMACH - 80085006  
 LIVER - 10200004  
 LUNG - 39607008  
 LYMPH NODE - 59441001  
 MEDIASTINUM - 72410000  
 NECK, LEFT SIDE - 170583000  
 NECK, RIGHT SIDE - 170303002  
 PANCREAS - 15776009  
 PAROTID GLAND - 45289007  
 PELVIS - 12921003  
 PERIRENAL TISSUE - 47145004  
 PERITONEAL CYTOLOGIC MATERIAL - 110944006  
 PERITONEUM - 15425007  
 PLEURAL CYTOLOGIC MATERIAL - 110913002  
 PROSTATE - 41216001  
 PYLORIC ANTRUM - 66051006  
 PYLORUS - 280119005  
 RECTOSIGMOID JUNCTION - 49832006  
 RECTUM - 34402009  
 RENAL PELVIS - 25990002  
 RETROPERITONEUM - 82849001  
 RIGHT COLIC FLEXURE - 48338005  
 RIGHT TESTIS - 15598003  
 SALIVARY GLAND - 385294005  
 SIGMOID COLON - 60184004  
 SKIN - 39937001

SMALL INTESTINE - 30315005  
 SOFT TISSUES - 87784001  
 SPLEEN - 78961009  
 STOMACH - 69695003  
 SUBLINGUAL GLAND - 88481005  
 SUBPHRENIC FOSSA - 243974009  
 SYNOVIAL CYTOLOGIC MATERIAL - 110895009  
 TAIL OF PANCREAS - 73239005  
 THYROID GLAND - 69748006  
 TISSUE - 85756007  
 TRANSVERSE COLON - 485005  
 URETER - 87953007  
 URINARY BLADDER - 89837001  
 URINE - 78014005  
 VAS DEFERENS - 57671007  
 VERTEBRA - 51282000

Follow the instructions outlined in the table below to perform the installation of Patch LR\*5.2\*469.

**Table 6: M Code Installation Instructions for Patch LR\*5.2\*469**

Step #	Description	✓
1	On the PackMan menu, use the <code>INSTALL/CHECK MESSAGE</code> option. This option loads the patch into a Transport Global on your system.	
2	From the Kernel Installation and Distribution System (KIDS) menu, select the <code>Installation</code> menu.	
3	From this menu, you may elect to use the following options:	
	<b>a</b> Verify Checksums in Transport Global: This option will allow you to ensure the integrity of the routines that are in the transport global.	
	<b>b</b> Print Transport Global: This option will allow you to view the components of the KIDS build.	
	<b>c</b> Compare Transport Global to Current System: This option will allow you to view all changes that will be made when this patch is installed. It compares all components of this patch (routines, Data Dictionaries (DD's), templates, etc.).	
	<b>d</b> Backup a Transport Global: This option will create a backup message of any routines exported with this patch. It will not backup any other changes such as DD's or templates.	
4	Under the <code>Installation</code> menu, select the package <code>LR*5.2*469</code> .	
5	When prompted 'Want KIDS to Rebuild Menu Trees Upon Completion of Install? <code>NO//'</code> , respond "NO".	

Step #	Description	✓
6	When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//', respond "NO".	
7	When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//', respond "NO".	
8	Enter the Device you want to print the Install message. You can queue the install by enter a "Q" at the device prompt. Enter a "^" to abort the install.	
9	When prompted 'Delay Install (Minutes): (0 - 60) ://', respond "0".	

## 4.13 Installation Verification Procedure

No specific verification procedures are required after the installations of patches LR\*5.2\*462 and LR\*5.2\*469. If desired, run Verify Checksums in Transport Global or use the following command: CHECK1^XTSUMBLD.

## 5 Back-Out Procedure

A Back-Out pertains to a return to the last known good operational state of the software and appropriate platform settings.

### 5.1 Back-Out Strategy

During installation, if the option to back up routines was run as directed, 'Backup a Transport Global', then routines have the ability to be restored from the "backup" MailMan message that was generated. However, the KIDS installation process does not perform a restore of other Vista components, such as data dictionary, cross-reference, and template changes, etc.

Prior to attempting a back-out of the software, contact the VA Help Desk at 1-855-673-4357 for support or assistance.

### 5.2 Back-Out Considerations

The LIM and the Chief of Pathology have the authority to order the back-out.

#### 5.2.1 Load Testing

Load testing was not performed for the patches.

### 5.3 Back-Out Criteria

Criteria for a back-out includes, but are not limited, to the following:

- Failed baseline testing.
- Non-recoverable software error.

## **5.4 Back-Out Risks**

No back-out risks have been determined at this time.

## **5.5 Authority for Back-Out**

The LIM and the Chief of Pathology have the authority to order the Back-out Procedure.

## **5.6 Back-Out Procedure**

Restore the backup message of any routines exported with the patches and reload the transport global.

## **5.7 Back-out Verification Procedure**

Verify the back-out procedure by comparing the transport global to the current system state.

# **6 Rollback Procedure**

Patches LR\*5.2\*462 and LR\*5.2\*469 and any installed dependent patches need to be taken out in reverse of the order in which these patches were installed; routines and data dictionary modifications and populated data must also be rolled back in reverse order.

Contact the VA Help Desk at 1-855-673-4357 for support or assistance regarding roll-back procedures.

## **6.1 Rollback Considerations**

No rollback considerations have been determined at this time.

## **6.2 Rollback Criteria**

The only criteria for a rollback that has been determined at this time is that the installation failed baseline testing.

## **6.3 Rollback Risks**

The only risk determined at this time is the possibility of downtime which would only effect the users of the Laboratory package.

## **6.4 Authority for Rollback**

The LIM and the Chief of Pathology have the authority to require the rollback and accept the risks.

## **6.5 Rollback Procedure**

The need for a rollback is highly unlikely, however if it is required, contact the product development team for assistance if needed. The rollback procedure may require Lab downtime and a reinstall of any previous KIDS versions.

## **6.6 Rollback Verification Procedure**

No specific rollback verification procedures are required after the installations of patches LR\*5.2\*462 and LR\*5.2\*469.

# Appendix A: LR\*5.2\*462 Captured Installation

The following is an installation example.

**Figure 10: LR\*5.2\*462 Captured Installation Example**

```
FRER1PSVR:FRE 21e1>D CHECK1^XTSUMBLD
```

New CheckSum CHECK1^XTSUMBLD:

This option determines the current checksum of selected routine(s).

The Checksum of the routine is determined as follows:

1. Any comment line with a single semi-colon is presumed to be followed by comments and only the line tag will be included.
2. Line 2 will be excluded from the count.
3. The total value of the routine is determined (excluding exceptions noted above) by multiplying the ASCII value of each character by its position on the line and position of the line in the routine being checked.

Select one of the following:

P	Package
B	Build

Build from: Package

All Routines? No => No

Routine: LR7OB3

Routine: LR7OF0

Routine: LR7OFAO

Routine: LRAP

Routine: LRAPDA

Routine: LRAPLG

Routine: LRAPLG1

Routine: LRAPRES1

```

Routine: LROS
Routine: LRUDEL
Routine: LRUWLF
Routine:
11 routines

LR7OB3      value = 24033857
LR7OF0      value = 40209802
LR7OFAO     value = 37973467
LRAP        value = 24889097
LRAPDA      value = 92935176
LRAPLG      value = 7969847
LRAPLG1     value = 169348471
LRAPRES1    value = 66194684
LROS        value = 33110400
LRUDEL      value = 25312700
LRUWLF      value = 10175468
done
FRER1PSVR:FRE 21e1>D ^XUP

Setting up programmer environment
This is a PRODUCTION account.

Terminal Type set to: C-VT220

You have 1774 new messages.
Select OPTION NAME: XPD MAIN      Kernel Installation & Distribution System

    Edits and Distribution ...
    Utilities ...
    Installation ...
    Patch Monitor Main Menu ...

Select Kernel Installation & Distribution System Option: Installation

    1      Load a Distribution
    2      Verify Checksums in Transport Global
    3      Print Transport Global

```



```

4      Compare Transport Global to Current System
5      Backup a Transport Global
6      Install Package(s)
      Restart Install of Package(s)
      Unload a Distribution

Select Installation Option: 2  Verify Checksums in Transport Global
Select INSTALL NAME:      LR*5.2*462      Loaded from Distribution      11/16/16@07:4
9:59

      => LR*5.2*462

This Distribution was loaded on Nov 16, 2016@07:49:59 with header of
      LR*5.2*462
It consisted of the following Install(s):
      LR*5.2*462

Want each Routine Listed with Checksums: Yes//      YES
DEVICE: HOME//      HOME      (CRT)

PACKAGE: LR*5.2*462      Nov 16, 2016 8:53 am      PAGE 1
-----
LR462A      Calculated      5945552
LR7OAPKM      Calculated      69415435
LR7OAPMG      Calculated      59502319
LR7OB3      Calculated      27184019
LR7OF0      Calculated      41838489
LR7OFAO      Calculated      54432906
LRAP      Calculated      25886264
LRAPDA      Calculated      95028265
LRAPKLG      Calculated      46698322
LRAPKLG1      Calculated      8921740
LRAPKOE      Calculated      164456408
LRAPKOE1      Calculated      1258642
LRAPKOEU      Calculated      7473474
LRAPLG      Calculated      12374267
LRAPLG1      Calculated      181070435
LRAPRES1      Calculated      70400068
Enter RETURN to continue or '^' to exit:

```

-----

LROS	Calculated	42665920
LRUDEL	Calculated	28725080
LRUWLF	Calculated	10789696

19 Routines checked, 0 failed.

- 1 Load a Distribution
  - 2 Verify Checksums in Transport Global
  - 3 Print Transport Global
  - 4 Compare Transport Global to Current System
  - 5 Backup a Transport Global
  - 6 Install Package(s)
- Restart Install of Package(s)
- Unload a Distribution

Select Installation Option: 5 Backup a Transport Global

Select INSTALL NAME: LR\*5.2\*462 Loaded from Distribution 11/16/16@07:49:59

=> LR\*5.2\*462

This Distribution was loaded on Nov 16, 2016@07:49:59 with header of

LR\*5.2\*462

It consisted of the following Install(s):

LR\*5.2\*462

Subject: Backup of LR\*5.2\*462 install on Nov 16, 2016

Replace

Loading Routines for LR\*5.2\*462.

Routine LR7OAPKM is not on the disk..

Routine LR7OAPMG is not on the disk.....

Routine LRAPKLG is not on the disk..

Routine LRAPKLG1 is not on the disk..

Routine LRAPKOE is not on the disk..

Routine LRAPKOE1 is not on the disk..

Routine LRAPKOEU is not on the disk.....

Send mail to: LASTNAME,FIRSTNAME

- 1 Load a Distribution
- 2 Verify Checksums in Transport Global
- 3 Print Transport Global
- 4 Compare Transport Global to Current System
- 5 Backup a Transport Global
- 6 Install Package(s)  
Restart Install of Package(s)  
Unload a Distribution

Select Installation Option: 6 Install Package(s)

Select INSTALL NAME: LR\*5.2\*462 Loaded from Distribution 11/16/16@07:49:59

=> LR\*5.2\*462

This Distribution was loaded on Nov 16, 2016@07:49:59 with header of

LR\*5.2\*462

It consisted of the following Install(s):

LR\*5.2\*462

Checking Install for Package LR\*5.2\*462

Will first run the Environment Check Routine, LR462A

--- Environment Check is Ok ---

Sending install loaded alert to mail group G.LMI

Install Questions for LR\*5.2\*462

Incoming Files:

60 LABORATORY TEST

Note: You already have the 'LABORATORY TEST' File.

69 LAB ORDER ENTRY

Note: You already have the 'LAB ORDER ENTRY' File.

69.71 LR CPRS PARAMETERS (including data)

```

69.9      LABORATORY SITE  (Partial Definition)
Note:  You already have the 'LABORATORY SITE' File.

Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//

Want KIDS to INHIBIT LOGONs during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME//      HOME  (CRT)

-----

Install Started for LR*5.2*462 :
      Nov 16, 2016@08:53:52

Build Distribution Date: Oct 06, 2016

Installing Routines:
      Nov 16, 2016@08:53:53

Running Pre-Install Routine: PRE^LR462A

      *** Preinstall completed ***

Installing Data Dictionaries:
      Nov 16, 2016@08:53:53

Installing Data:
      Nov 16, 2016@08:53:53

Installing PACKAGE COMPONENTS:

Installing PRINT TEMPLATE

```

Installing INPUT TEMPLATE

Installing OPTION

Nov 16, 2016@08:53:53

Running Post-Install Routine: POST^LR462A

Sending install completion alert to mail group G.LMI

Updating Routine file...

Updating KIDS files...

LR\*5.2\*462 Installed.

Nov 16, 2016@08:53:54

LR\*5.2\*462

Install Completed

## Appendix B: LR\*5.2\*469 Captured Installation

The following is an installation example.

**Figure 11: LR\*5.2\*469 Captured Installation Example**

```
Setting up programmer environment
This is a PRODUCTION account.

Terminal Type set to: C-VT220

You have 5565 new messages.
Select OPTION NAME: XPD MAIN      Kernel Installation & Distribution System

    Edits and Distribution ...
    Utilities ...
    Installation ...
    Patch Monitor Main Menu ...

Select Kernel Installation & Distribution System Option: Installation

    1      Load a Distribution
    2      Verify Checksums in Transport Global
    3      Print Transport Global
    4      Compare Transport Global to Current System
    5      Backup a Transport Global
    6      Install Package(s)
           Restart Install of Package(s)
           Unload a Distribution
           Check Local Copies of Incoming Routines
           Find Active Options
           Only Lookup in Install File

Select Installation Option: 2  Verify Checksums in Transport Global
Select INSTALL NAME:      LR*5.2*469      11/22/16@15:12:05
=> LR*5.2*469

This Distribution was loaded on Nov 22, 2016@15:12:05 with header of
    LR*5.2*469 TEST v10
    It consisted of the following Install(s):
```

LR\*5.2\*469

Want each Routine Listed with Checksums: Yes// YES

DEVICE: HOME// HOME (CRT)

PACKAGE: LR\*5.2\*469 Nov 22, 2016 3:37 pm

PAGE 1

-----  
LR469A Calculated 59641039

1 Routine checked, 0 failed.

- 1 Load a Distribution
  - 2 Verify Checksums in Transport Global
  - 3 Print Transport Global
  - 4 Compare Transport Global to Current System
  - 5 Backup a Transport Global
  - 6 Install Package(s)
- Restart Install of Package(s)
- Unload a Distribution
- Check Local Copies of Incoming Routines
- Find Active Options
- Only Lookup in Install File

Select Installation Option: 5 Backup a Transport Global

Select INSTALL NAME: LR\*5.2\*469 11/22/16@15:12:05

=> LR\*5.2\*469

This Distribution was loaded on Nov 22, 2016@15:12:05 with header of

LR\*5.2\*469

It consisted of the following Install(s):

LR\*5.2\*469

Subject: Backup of LR\*5.2\*469 install on Nov 22, 2016

Replace

Loading Routines for LR\*5.2\*469.

Send mail to: LASTNAME,FIRSTNAME

- 1 Load a Distribution
- 2 Verify Checksums in Transport Global

```

3      Print Transport Global
4      Compare Transport Global to Current System
5      Backup a Transport Global
6      Install Package(s)
      Restart Install of Package(s)
      Unload a Distribution
      Check Local Copies of Incoming Routines
      Find Active Options
      Only Lookup in Install File

Select Installation Option: 6  Install Package(s)

DON'T FORGET to use 'Check Local Copies of Incoming Routines' Option FIRST!

Select INSTALL NAME:      LR*5.2*469      11/22/16@15:12:05
      => LR*5.2*469

This Distribution was loaded on Nov 22, 2016@15:12:05 with header of
      LR*5.2*469
      It consisted of the following Install(s):
      LR*5.2*469
Checking Install for Package LR*5.2*469
Will first run the Environment Check Routine, LR469A

      Found ORDERABLE ITEM 'BONE MARROW'
      Found ORDERABLE ITEM 'BRONCHIAL BIOPSY'
      Found ORDERABLE ITEM 'BRONCHIAL CYTOLOGY'
      Found ORDERABLE ITEM 'DERMATOLOGY'
      Found ORDERABLE ITEM 'FINE NEEDLE ASPIRATE'
      Found ORDERABLE ITEM 'GASTROINTESTINAL ENDOSCOPY'
      Found ORDERABLE ITEM 'GENERAL FLUID'
      Found ORDERABLE ITEM 'GYNECOLOGY (PAP SMEAR) '
      Found ORDERABLE ITEM 'TISSUE EXAM'
      Found ORDERABLE ITEM 'UROLOGY,BLADDER/URETER'
      Found ORDERABLE ITEM 'UROLOGY,PROSTATE'
      Found ORDERABLE ITEM 'RENAL BIOPSY'
      Found ORDERABLE ITEM 'URINE '

      *** 'ABDOMEN'

```



```
Found SPECIMEN for SNOMED CT '113345001'

*** 'ADRENAL GLAND'
Found SPECIMEN for SNOMED CT '23451007'

*** 'ANAL CANAL'
Found SPECIMEN for SNOMED CT '34381000'

*** 'APPENDIX'
Found SPECIMEN for SNOMED CT '66754008'

*** 'ASCENDING COLON'
Found SPECIMEN for SNOMED CT '9040008'

*** 'BILE DUCT CYTOLOGIC MATERIAL'
Found SPECIMEN for SNOMED CT '110928002'

*** 'BILE DUCT MUCOUS MEMBRANE'
Found SPECIMEN for SNOMED CT '7035006'

*** 'BILIARY TRACT'
Found SPECIMEN for SNOMED CT '34707002'

*** 'BODY OF PANCREAS'
Found SPECIMEN for SNOMED CT '40133006'

*** 'BONE'
Found SPECIMEN for SNOMED CT '3138006'

*** 'BONE MARROW'
Found SPECIMEN for SNOMED CT '14016003'

*** 'BREAST'
Found SPECIMEN for SNOMED CT '76752008'

*** 'BRONCHIAL CYTOLOGIC MATERIAL'
Found SPECIMEN for SNOMED CT '110912007'

*** 'BRONCHUS'
```

```
Found SPECIMEN for SNOMED CT '955009'

*** 'CARDIAC INCISURE OF STOMACH'
    Found SPECIMEN for SNOMED CT '5459006'

*** 'CARDIAC OSTIUM OF STOMACH'
    Found SPECIMEN for SNOMED CT '63853002'

*** 'CARDIO-ESOPHAGEAL JUNCTION'
    Found SPECIMEN for SNOMED CT '25271004'

*** 'CECUM'
    Found SPECIMEN for SNOMED CT '32713005'

*** 'CEREBROSPINAL FLUID CYTOLOGIC MATERIAL'
    Found SPECIMEN for SNOMED CT '110969006'

*** 'CERVICAL CYTOLOGIC MATERIAL'
    Found SPECIMEN for SNOMED CT '110949001'

*** 'CHORIONIC VILLI'
    Found SPECIMEN for SNOMED CT '2049008'

*** 'COLON'
    Found SPECIMEN for SNOMED CT '71854001'

*** 'DESCENDING COLON'
    Found SPECIMEN for SNOMED CT '32622004'

*** 'DUODENUM'
    Found SPECIMEN for SNOMED CT '38848004'

*** 'ESOPHAGUS'
    Found SPECIMEN for SNOMED CT '32849002'

*** 'ESOPHAGUS, LOWER THIRD'
    Found SPECIMEN for SNOMED CT '67173009'

*** 'ESOPHAGUS, MIDDLE THIRD'
```

```
Found SPECIMEN for SNOMED CT '19000002'

*** 'ESOPHAGUS, UPPER THIRD'
Found SPECIMEN for SNOMED CT '54738009'

*** 'GASTRIC FUNDUS'
Found SPECIMEN for SNOMED CT '414003'

*** 'GASTRIC JUICE'
Found SPECIMEN for SNOMED CT '31773000'

*** 'GREATER CURVATURE OF STOMACH'
Found SPECIMEN for SNOMED CT '89382009'

*** 'HEAD OF PANCREAS'
Found SPECIMEN for SNOMED CT '64163001'

*** 'ILEUM'
Found SPECIMEN for SNOMED CT '34516001'

*** 'JEJUNUM'
Found SPECIMEN for SNOMED CT '21306003'

*** 'KIDNEY'
Found SPECIMEN for SNOMED CT '64033007'

*** 'LEFT COLIC FLEXURE'
Found SPECIMEN for SNOMED CT '72592005'

*** 'LEFT TESTIS'
Found SPECIMEN for SNOMED CT '63239009'

*** 'LESSER CURVATURE OF STOMACH'
Found SPECIMEN for SNOMED CT '80085006'

*** 'LIVER'
Found SPECIMEN for SNOMED CT '10200004'

*** 'LUNG'
```

```
Found SPECIMEN for SNOMED CT '39607008'

*** 'LYMPH NODE'
    Found SPECIMEN for SNOMED CT '59441001'

*** 'MEDIASTINUM'
    Found SPECIMEN for SNOMED CT '72410000'

*** 'NECK, LEFT SIDE'
    Found SPECIMEN for SNOMED CT '170583000'

*** 'NECK, RIGHT SIDE'
    Found SPECIMEN for SNOMED CT '170303002'

*** 'PANCREAS'
    Found SPECIMEN for SNOMED CT '15776009'

*** 'PAROTID GLAND'
    Found SPECIMEN for SNOMED CT '45289007'

*** 'PELVIS'
    Found SPECIMEN for SNOMED CT '12921003'

*** 'PERICARDIAL CYTOLOGIC MATERIAL'
    Found SPECIMEN for SNOMED CT '110919003'

*** 'PERIRENAL TISSUE'
    Found SPECIMEN for SNOMED CT '47145004'

*** 'PERITONEAL CYTOLOGIC MATERIAL'
    Found SPECIMEN for SNOMED CT '110944006'

*** 'PERITONEUM'
    Found SPECIMEN for SNOMED CT '15425007'

*** 'PLEURAL CYTOLOGIC MATERIAL'
    Found SPECIMEN for SNOMED CT '110913002'

*** 'PROSTATE'
```

```
Found SPECIMEN for SNOMED CT '41216001'

*** 'PYLORIC ANTRUM'
    Found SPECIMEN for SNOMED CT '66051006'

*** 'PYLORUS'
    Found SPECIMEN for SNOMED CT '280119005'

*** 'RECTOSIGMOID JUNCTION'
    Found SPECIMEN for SNOMED CT '49832006'

*** 'RECTUM'
    Found SPECIMEN for SNOMED CT '34402009'

*** 'RENAL PELVIS'
    Found SPECIMEN for SNOMED CT '25990002'

*** 'RETROPERITONEUM'
    Found SPECIMEN for SNOMED CT '82849001'

*** 'RIGHT COLIC FLEXURE'
    Found SPECIMEN for SNOMED CT '48338005'

*** 'RIGHT TESTIS'
    Found SPECIMEN for SNOMED CT '15598003'

*** 'SALIVARY GLAND'
    Found SPECIMEN for SNOMED CT '385294005'

*** 'SIGMOID COLON'
    Found SPECIMEN for SNOMED CT '60184004'

*** 'SKIN'
    Found SPECIMEN for SNOMED CT '39937001'

*** 'SMALL INTESTINE'
    Found SPECIMEN for SNOMED CT '30315005'

*** 'SOFT TISSUES'
```

```
Found SPECIMEN for SNOMED CT '87784001' and name 'SOFT TISSUES'

*** 'SPLEEN'
    Found SPECIMEN for SNOMED CT '78961009'

*** 'STOMACH'
    Found SPECIMEN for SNOMED CT '69695003'

*** 'SUBLINGUAL GLAND'
    Found SPECIMEN for SNOMED CT '88481005'

*** 'SUBMAXILLARY GLAND'
    Found SPECIMEN for SNOMED CT '54019009'

*** 'SUBPHRENIC FOSSA'
    Found SPECIMEN for SNOMED CT '243974009'

*** 'SYNOVIAL CYTOLOGIC MATERIAL'
    Found SPECIMEN for SNOMED CT '110895009'

*** 'TAIL OF PANCREAS'
    Found SPECIMEN for SNOMED CT '73239005'

*** 'THYROID GLAND'
    Found SPECIMEN for SNOMED CT '69748006'

*** 'TISSUE'
    Found SPECIMEN for SNOMED CT '85756007'

*** 'TRANSVERSE COLON'
    Found SPECIMEN for SNOMED CT '485005'

*** 'URETER'
    Found SPECIMEN for SNOMED CT '87953007'

*** 'URINARY BLADDER'
    Found SPECIMEN for SNOMED CT '89837001'

*** 'URINE'
```

```
Found SPECIMEN for SNOMED CT '78014005' and name 'URINE'

*** 'VAS DEFERENS'
Found SPECIMEN for SNOMED CT '57671007'

*** 'VERTEBRA'
Found SPECIMEN for SNOMED CT '51282000' and name 'VERTEBRA'

Install Questions for LR*5.2*469

Incoming Files:

69.73      LR AP DIALOG CONFIG  (including data)

Want KIDS to INHIBIT LOGONS during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME// QUEUE TO PRINT ON
DEVICE: HOME// P-M

1 P-MESSAGE-HFS  OUTPUT TO MAIL
2 P-MSG  MAILMAN
Choose 1-2> 1  P-MESSAGE-HFS  OUTPUT TO MAIL

Subject: LR*5.2*469T10 INSTALL LIVE - MAC

Select one of the following:

M          Me
P          Postmaster

From whom: Me//
Send mail to: LASTNAME,FIRSTNAME// FIRSTAME.LASTNAME@VA.GOV  FORUM.VA.GOV via
FORUM.VA.GOV
```

Request Start Time: Nov 22, 2016@15:40// 11/23@0600 (NOV 23, 2016@06:00:00)

Install Queued!

- 1      Load a Distribution
  - 2      Verify Checksums in Transport Global
  - 3      Print Transport Global
  - 4      Compare Transport Global to Current System
  - 5      Backup a Transport Global
  - 6      Install Package(s)
- Restart Install of Package(s)
- Unload a Distribution
- Check Local Copies of Incoming Routines
- Find Active Options
- Only Lookup in Install File